# GIMA

REF 28202-28232-28233-28234

**ASPIRATORE TOBI HOSPITAL**  
**ASPIRATOR TOBI HOSPITAL**  
**ASPIRATEUR TOBI HOSPITAL**  
**ABSAUGER TOBI HOSPITAL**  
**ASPIRADOR TOBI HOSPITAL**

---

**MANUALE D'USO**  
**INSTRUCTION MANUAL**  
**MODE D'EMPLOI**  
**HANDBUCH**  
**MANUAL DE ISTRUCCIONES**

**ATTENZIONE:** Gli operatori devono leggere e capire completamente questo manuale prima di utilizzare il prodotto.  
**ATTENTION:** The operators must carefully read and completely understand the present manual before using the product.  
**AVIS:** Les opérateurs doivent lire et bien comprendre ce manuel avant d'utiliser le produit.  
**ACHTUNG:** Die Bediener müssen vorher dieses Handbuch gelesen und verstanden haben, bevor sie das Produkt benutzen.  
**ATENCIÓN:** Los operadores tienen que leer y entender completamente este manual antes de utilizar el producto.

CE 0476

M28202-M-Rev.6.11.12

**GIMA S.p.A.**  
Via Marconi, 1 – 20060 Gessate (MI) Italia  
**ITALIA:** Tel. 02 953854.1 – Fax: 02 95381167  
E-mail: [gima@gimaitaly.com](mailto:gima@gimaitaly.com) –  
[www.gimaitaly.com](http://www.gimaitaly.com)

**INTERNATIONAL:** Tel. +39 02 953854209/221/225  
Fax. +39 02 95380056  
E-mail: [export@gimaitaly.com](mailto:export@gimaitaly.com) – [www.gimaitaly.com](http://www.gimaitaly.com)

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**TOBI HOSPITAL** is a surgical aspirator working with 230V ~ / 50 Hz network electricity, to be used for the aspiration of the body liquids (mucus or catarrh and blood). The device is equipped with 4 antistatic wheels, two with brakes, and movement handle. Easily movable and designed for virtually continuous use. Thanks to this characteristics and to the rating that it has, this product is particularly suitable for hospital use, on the tracheotomized patients and surgical applications. Made of highly heat resistant, electrically insulated plastic material in conformity with the latest European safety standard, the product is supplied with two polycarbonate autocavable jars with overflow valve and equipped with aspiration regulator and vacuum indicator located on the front panel. Available under request with different versions for applications and use (version with remote control, version with 4 litre jars, version with flow direction regulator).

## GENERAL WARNING

**READ INSTRUCTION MANUAL CAREFULLY BEFORE USE**

**ONLY HIGHLY QUALIFIED STAFF USE RESERVED**

**THE INSTRUMENT MUST NOT BE DISASSEMBLED. FOR A TECHNICAL SERVICE ALWAYS CONTACT GIMA SpA**

**KEEP OFF THE REACH OF CHILDREN OR NOT CAPABLE PEOPLE WITHOUT SUPERVISION**

**FULL CONTAINERS MUST BE HANDLED WITH GREAT CARE DURING TRANSFER TO THE DISPOSAL AREAS, FOLLOWING THE LOCAL PROCEDURES AND REGULATIONS**

## IMPORTANT SAFETY RULES

1. On opening the packaging, check the integrity of the appliance, paying particular attention to the presence of damage to the plastic parts, which may make access possible to internal live parts and also to breakage and/or peeling of the power supply cable. **In these cases don't connect the plug to the electric socket. Carry out these controls before each use.**
2. before connecting the appliance always check that the electric data indicated on the data label and the type of plug used, correspond to those of the mains electricity to which it's to be connected.
3. If the plug supplied with the appliance is incompatible with the mains electricity socket, contact qualified staff for replacement of the plug with a suitable type. The use of simple or multiple and/or extension adapters is not generally recommended. Whenever their use is indispensable, use those in compliance with safety regulations, however paying attention not to exceed the maximum power supply limits, which are indicated on the adapters and extensions.
4. Respect the safety regulations indicated for electrical appliances and particularly:
  - Use original components and accessories provided by the manufacturer to guarantee the highest efficiency and safety of the device.
  - The device can be used only with the bacteriological filter.
  - Never immerse the appliance into water.
  - Avoid touching the aspirator with wet hands and always prevent the equipment from getting in touch with liquids. Never leave the equipment near water or immerse it into a liquid. Should the equipment fall into water, detach its power cable from the socket before touching it.
  - None of the electrical and/or mechanical parts of the machine is designed to be repaired by the client and/or by its user. Do not open the aspirator or disassembly its electrical and/or mechanical parts. Always report to GIMA SpA technical support.
  - Using the equipment in environmental conditions other than those indicated in this manual may seriously endanger its safety and technical parameters.
  - Position the appliance on flat stable surfaces.
  - Position the device in a way that the air inlets on the back aren't obstructed.
  - Never use the device in environments which have anaesthetic mixtures inflammable with air, oxygen or nitric oxide.
  - Don't touch the device with wet hands and always prevent the appliance coming into contact with liquids.
  - Keep off the reach of children or not capable people without supervision.
  - Don't leave the appliance connected to the power supply socket when not in use.
  - Don't pull the power supply cable to disconnect the plug remove the plug from the mains socket correctly.
  - Preserve and use the medical device in environments protected from atmospheric factors and at a distance from heat sources.
  - Don't use the device thoracic drainage.
5. For repairs, exclusively contact GIMA SpA technical service and request the use of original spare parts. Failure to comply with the above can jeopardise the safety of the device.
6. **This medical device must be destined exclusively for the use for which it has been designed and described in this manual.** Any different use must be considered incorrect and therefore dangerous; the manufacturer cannot be considered liable for damage caused by improper, incorrect and / or unreasonable use or if the appliance is used in electrical plants that are not in compliance with the regulations in force.
7. Particular precautions must be made concerning electromagnetic compatibility. The medical device must be installed and used according to information supplied with the accompanying documents.
8. Instrument and accessories discharging must be done following current law regulations in every country of use.
9. None of electric or mechanical parts have been designed to be repaired by customers or end-users. Don't open the device, do not mishandle the electric / mechanical parts. Always contact GIMA SpA technical assistance.
10. Using the device in environmental conditions different than those indicated in this manual may harm seriously the safety and the technical characteristics of the same.

The manufacturer cannot be held liable for accidental or indirect damages should the device be modified, repaired without authorization or should any of its component be damaged due to accident or misuse.

Any minimal modification / repair on the device voids the warranty and does not guarantee the compliance with the technical requirements provided by the MDD 93/42/EEC Directive and its normatives.

### TECHNICAL CHARACTERISTICS

TYPOLOGY (MDD 93/42/EEC)	Class IIa Medical Decice
MODEL	<b>TOBI HOSPITAL</b>
UNI EN ISO 10079-1	HIGH VACUUM / HIGH FLOW
POWER FEEDING	230V~ / 50Hz
POWER CONSUMPTION	300 VA
FUSE	F 1 x 4A L 250V
MAXIMUM SUCTION PRESSURE (without jar)	-90 kPa (-0.90 Bar)
MAXIMUM SUCTION FLOW (without jar)	90 l/min
WEIGHT	20 Kg
SIZE	650 x 550 x 450 mm
DUTY CYCLE (to 35°C and 110% operating voltage)	Non-Stop Operated
SICILICONE TUBE SIZE	Ø 8 x 14 mm
ACCURANCY OF VACUUM INDICATOR	± 5%
WORKING CONDITION	Room temperature: 5 ÷ 35°C Room humidity percentage: 30 ÷ 75% RH Altitude: 0 ÷ 2000m s.l.m.
CONSERVATION CONDITION AND TRASPORT	Room temperature: -40 ÷ 70°C Room humidity percentage: 10 ÷ 100% RH

### SYMBOLS

	Class II isolation equipment
	CE marking in conformity with EC directive 93/42/EEC and subsequent changes. Manufactured by: GIMA S.P.A. - Via Marconi, 1 - 20060 Gessate (MI) - Italy
	Warning, consult the instruction manual
	To Preserve in place coolness and dry land
	Conservation temperature: -40 ÷ 70°C
	Type B equipment
	Fuse
~	Alternate Current
Hz	Mains Frequency
I	ON
O	OFF
II	Remote Control

Please note technical specifications may vary upon the manufacturer's discretion!

### CLEANING THE MAIN UNIT

To clean the device external parts always use a cotton cloth dampened with detergent.  
Don't use abrasive or solvent detergents.

**PARTICULAR CARE SHOULD BE TAKEN TO ENSURE THAT THE INTERNAL PARTS OF THE EQUIPMENT DO NOT GET IN TOUCH WITH LIQUIDS. NEVER CLEAN THE EQUIPMENT WITH WATER.**

During all clearing operations use protection gloves and apron (if need be, also wear a face mask and glasses) to avoid getting in contact with contaminating substances (after each utilization cycle of the machine).

Guidance and manufacturer's declaration – Electromagnetic Emissions		
The surgical aspirator <b>TOBI HOSPITAL</b> is intended for use in the electromagnetic environment specified below. The customers or the user of the surgical aspirator <b>TOBI HOSPITAL</b> should assure that it's used in such an environment.		
Emissions Test	Compliance	Electromagnetic environment - guidance
Irradiated / Conducted emissions CISPR11	Group 1	The surgical aspirator TOBI HOSPITAL only used RF energy only for its internal functioning. Therefore its RF emissions are very low and are not cause interference in proximity of any Electronic appliances.
Irradiated / Conducted emissions CISPR11	Class [B]	The surgical aspirator TOBI HOSPITAL can be used in all environments, including domestic and those connected directly to the public mains distribution that supplies power to environments used for domestic scopes.
Harmonic emissions IEC/EN 61000-3-2	Class [A]	
Voltage fluctuations / flicker emissions IEC/EN 61000-3-3	Complies	

Guidance and manufacturer's declaration – Immunity Emissions			
The surgical aspirator <b>TOBI HOSPITAL</b> is intended for use in the electromagnetic environment specified below. The customers or the user of the surgical aspirator <b>TOBI HOSPITAL</b> should assure that it's used in such an environment.			
Immunity Test	Level indicated by the IEC 60601-1-2	Compliance Level	Electromagnetic environments - guidance
Electrostatic discharge (ESD) IEC/EN 61000-4-2	± 6kV on contact ± 8kV in air	The device doesn't change its state	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC/EN 61000-4-4	± 2kV power supply lines  ± 1kV for input / output lines	The device doesn't change its state	Mains power quality should be that of a typical commercial environment or hospital.
Surge IEC/EN 61000-4-5	± 1kV differential mode	The device doesn't change its state	Mains power quality should be that of a typical commercial environment or hospital.
Loss of voltage, brief voltage interruptions and variations IEC/EN 61000-4-11	5%U <sub>T</sub> (>95% dip U <sub>T</sub> ) for 0.5 cycle 40%U <sub>T</sub> (>60% dip U <sub>T</sub> ) for 5 cycle 70%U <sub>T</sub> (>30% dip U <sub>T</sub> ) for 25 cycle <5%U <sub>T</sub> (>95% dip U <sub>T</sub> ) for 5 sec	-	Mains power quality should be that of a typical commercial environment or hospital. If the user of the surgical aspirator TOBI HOSPITAL request that the appliance operates continuously, the use of a continuity unit is recommended.
Magnetic field IEC/EN 61000-4-8	3A/m	The device doesn't change its state	The power frequency magnetic field should be measured in the intended installation location to assure that it's sufficiently low.
Note U <sub>T</sub> is the value of the power supply voltage			

Guidance and manufacturer's declaration – Immunity Emissions			
The surgical aspirator <b>TOBI HOSPITAL</b> is intended for use in the electromagnetic environment specified below. The customers or the user of the surgical aspirator <b>TOBI HOSPITAL</b> should assure that it's used in such an environment.			
Immunity Test	Level indicated by the IEC 60601-1-2	Compliance level	Electromagnetic environments - guidance
Conducted Immunity IEC / EN 61000-4-6	3Vrms 150kHz to 80MHz (for non life-supporting devices)	V <sub>1</sub> = 3 V rms	<p>The portable and mobile RF communication devices, including cables, must not be used closer to the TOBI HOSPITAL device, than the separation distance calculated by the equation applicable to the transmitter frequency. Recommended separation distance</p> $d = [3.5 / V_1] \sqrt{P}$ $d = [12 / E_1] \sqrt{P} \text{ from 80 MHz to 800MHz}$ $d = [23 / E_1] \sqrt{P} \text{ from 800 MHz to 2.5 GHz}$ <p>Where P is the maximum nominal output voltage of the transmitter in Watt (W) depending on the manufacturer of the transmitter and the recommended separation distance in metres (m). The intensity of the field from the fixed RF transmitters, as determined by an electro-magnetic study of the site<sup>a)</sup>, could be lower than the level of conformity of each frequency interval<sup>b)</sup>. It is possible to check for interference in proximity to devices identified by the following symbol:</p> 
Radiated Immunity IEC / EN 61000-4-3	3V/m 80MHz to 2.5GHz (for non life-supporting devices)	E <sub>1</sub> = 3 V / m	
Note 1: At 80 MHz and 800 MHz the interval with the highest frequency is applied Note 2: These guide lines may not be applicable in all situations. The electro-magnetic propagation is influenced by the absorption and by			

reflection from buildings, objects and people.

a) The field intensity for fixed transmitters such as the base stations for radiotelephones (mobile and cordless) and terrestrial mobile radio, amateur radio devices, radio AM and FM transmitters and TV transmitters can not be theoretically and accurately foreseen.

To establish an electro-magnetic environment generated by fixed RF transmitters, an electro-magnetic study of the site should be considered. If the field intensity measured in the place where the device will be used surpasses the above mentioned applicable level of conformity, the normal functioning of the device should be monitored. If abnormal performance arises, additional measures such as changing the device's direction or positioning may be necessary.

b) The field intensity on an interval frequency of 150 kHz to 80 MHz should be less than 3 V/m.

#### Recommended separation distance between portable and mobile radio-communication devices and the monitor

The TOBI HOSPITAL surgical aspirator is intended to operate in an electro-magnetic environment where RF irradiated interferences are under control. The client or operator of the TOBI HOSPITAL device can help prevent electro-magnetic interference by keeping a minimum distance between the portable and mobile RF communication devices (transmitters) and the TOBI HOSPITAL device, as recommended below, in relation to the radio-communication maximum output power.

Maximum nominal output power of the Transmitter W	Separation distance from the frequency transmitter (m)		
	150 kHz to 80 MHz $d = [3.5 / V_1] \sqrt{P}$	80 MHz to 800 MHz $d = [12/E_1] \sqrt{P}$	800 MHz to 2.5 GHz $d = [23/E_1] \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters with a maximum nominal output power not shown above, the recommended separation distance in metres (m) can be calculated using the equation applicable to the transmitter frequency, where P is the maximum nominal output power of the transmitter in Watt (W) depending on the transmitter's manufacturer.

Note 1: At 80 MHz and 800 MHz the interval with the highest frequency is applied

Note 2: These guide lines may not be applicable in all situations. The electro-magnetic propagation is influenced by the absorption and by the reflection from buildings, objects and people.

This section contains information regarding the conformity of the compliance with the IEC 60601-1-2 Standard.

The TOBI HOSPITAL surgical aspirator is an electro-medical device that requires particular precautions regarding electro-magnetic compatibility and which must be installed and commissioned according to the electro-magnetic compatibility information supplied.

Mobile and portable RF communication appliances (mobile phones, transceivers, etc..) can affect the medical system.

The use of accessories, transducers and cables different to those specified, with the exception of transducers and cables sold by the appliance and system manufacturer as spare parts, can lead to an increase in emissions or in a decrease of the immunity of the device or system.

### ACCESSORIES SUPPLIED

DESCRIPTION
N°2 COMPLETE ASPIRATION JAR 2000cc
CONICAL FITTING
TUBES SET 8 mm x 14 mm
ANTIBACTERIAL AND HYDROFOBIC FILTER
FOOTSWITCH CONTROL (for versions equipped with footswitch control)

#### Replacing the antibacterial filter:

The filter is made of hydrophobic material that stops the passage of liquids into the same filter.

If you suspect the filter may have been contaminated and/or got wet or discoloured, always remove and replace the filter.

If the equipment is to be used on patients with unknown pathological conditions or should you evaluate the possibility of indirect contamination, remove and **replace the filter after each utilization**. The filter is not designed for decontamination, disassembly and/or sterilization. If you suspect the filter may have been contaminated and/or got wet or discoloured, always remove and replace the filter. If the equipment is to be used on patients whose pathologies are known and not implying any indirect contamination risks, we recommend to remove and replace the filter at the end of each work shift or else every month, even if the equipment has not been used. 4000cc complete tank versions are available on request.

Versions fitted with FLOVAC® 2000ml or 3000ml disposable collection systems (including a re-usable rigid polycarbonate container and a disposable Liner) are also available on request.

**WARNING:** Suction tubes for insertion in the human body purchased separately from the machine should comply with ISO 10993-1 standards on material biocompatibility.

## CLEANING ACCESSORIES AND INTERNAL PARTS

At the end of the application switch the equipment off and clean all its accessories as follows:

- Wear protection gloves and apron (if need be, also wear a face mask and glasses) to avoid getting in contact with contaminating substances;
- Disconnect the tank from the equipment removing any tubes connected to the container and paying particular attention to avoiding accidental contaminations.
- Empty and dispose of the flacon content complying with hospital regulations as well as with any provisions in force, including local regulations.
- Separate all the parts of the lid (float device and rings).

After disposing of disposable parts and disassembling the jar wash in running cold water and rinse thoroughly.

Then soak in warm water (temperature shall not exceed 60°C). Wash thoroughly and if necessary use a non-abrasive brush to remove incrustations. Rinse in running warm water and dry all parts with a soft cloth (non-abrasive).

The jar and the cover can be autoclaved by placing the parts into the autoclave and running one sterilization stem cycle at 121°C (1 bar relative pressure) making sure that the jar is positioned upsidedown. Mechanical resistance of the jar is guaranteed up to 30 cycles of sterilization and cleaning at the indicated conditions (EN ISO 10079-1). Beyond this limit the physical-mechanical characteristics of the plastic may decrease and replacement of the part is therefore recommended.

After sterilization and cooling at environment temperature of the parts make sure that these are not damaged.

Assemble the jar as follows:

- Place the overflow valve into its seat in the cover (under VACUUM connector)
- Insert floating valve keeping the o-ring towards the opening of the cage
- Place the o-ring into its seat around the cover
- After completing assembling operations always make sue that cover seals perfectly to avoid vacuum leakages or liquid exit

The aspiration tubes can be sterilized on autoclave using a sterilization cycle at 120°C.

The conical connector can be sterilized on autoclave using a sterilization cycle at 121°C.

The device is ready for a new employment now.

**DO NOT WASH, STERILIZE OR PUT IN AUTOCLAVE THE ANTIBACTERIAL FILTER**

### **Instruction for disposal Liner Flovac®:**

If the device is equipped with disposable collection systems FLOVAC ® carry out the disposal of the bag as follows:

Turn off the Vacuum and remove all the tubes connected to the Liner, giving particular attention to avoid accidental contamination. Fit the appropriate plugs to the "PATIENT" and "TANDEM" ports, pressing the home firmly, taking care to avoid accidental contamination. Turn the butterfly connector to OFF. Remove the liner bag from the rigid container and transfer it to the waste disposal area, ensuring that all the openings are sealed, keeping in mind the product is potentially infectious. This product must be disposed of in accordance with the current hospital regulations.

## RULES FOR RETURNING AND REPAIRING

**COMPLYING WITH THE NEW EUROPEAN RULES, GIMA INDICATES THE IMPORTANT POINTS TO PROTECT INSTRUMENT AND OPERATORS HYGIENE. THESE RULES MUST BE RESPECTED IN ORDER TO GUARANTEE HYGIENE AND SAFETY TO ALL THE PEOPLE OPERATING WITH THE INSTRUMENT TO OBTAIN QUALITY AND WELL BEING.**

GIMA warrants it's products for **24 months** after purchasing date.

In front of this warranty, GIMA will be obliged only to repair or substitute free of charge the products or parts of them that, after verification effected on our factory, or our authorized Service Center, by the Technical Service, results defective.

The product must be accompanied by a description of the defect.

The warranty, with exclusion of responsibility for direct and indirect damages, is thought limited to the solos defects of material or workmanship and it stops having effect when the device results however gotten off, tampered or sheltered out of the Factory or from the Authorized Service center. The commodity always travels to risk and danger of the buyer, without any responsibility of GIMA for damages caused by the transport or dismay from the vector. Every returned instrument will be hygienically checked before repairing.

If GIMA finds instrument not suitable for repairing due to clear signs of internal or external contamination, the same will be returned to customer with specification of NOT REPAIRED INSTRUMENT, accompanied by an explanation letter.

GIMA will decide if contamination is due to bad functioning or misuse. If contamination is due to bad functioning, GIMA will substitute the instrument, only if a SALE RECEIPT and STAMPED GUARANTEE accompany the same.

GIMA is not responsible for contaminated accessories, they will be substitute at customer's expenses.

For this reason it is **COMPULSORY** to carefully disinfect the external part of the instrument and accessories with a cloth soaked in methylated spirits or hypochlorite-based solutions. Put the instrument and accessories in a bag with indication of disinfecting. We also request to specify the kind of fault, in order to speed up repairing procedures.

To this end, please read the instructions carefully in order to avoid damaging the equipment through improper use.

Always specify the fault encountered so that GIMA can establish whether it falls into the category of the faults covered by the guarantee.

## MAINTENANCE

The **TOBI HOSPITAL** suction equipment does not need maintenance or lubrication.

It is necessary to check functioning and instrument before every use.

Unpack the instrument and **always check** integrity of plastic parts and feeding cable, they might have been damaged during previous use. Connect cable to electrical network and turn switch on.

Close the aspiration outlet with your finger and with suction regulator in maximum vacuum position check that the vacuum indicators reaches -90 kPa (-0.90 bar) maximum. Rotate the knob from right to left and check the aspiration regulating control.

The vacuum indicator should go down -40 kPa (-0.40 bar).

Verify that loud noises are not present, these can indicate wrong functioning.

A protection fuse (**F 1 x 4A L 250V**) reachable from exterior and it situated in the plug protects the instrument. For use replacing, always check the type and the range indicated.

Fault type	Cause	Solution
1. The suction unit doesn't work	Cable is damaged	Replace the cable
	External power source failure	Check the external power source
2. No aspiration	Jar Cap badly screwed down	Unscrewed the cap, then rescrew it correctly
3. No aspiration	Lid seal not in its seat	Unscrew the cap and insert the seal properly in its seat
4. The Vacuum power on the patient side is either very low or absent	a) Vacuum regulator set to minimum	a) Turn the vacuum regulator clockwise and check the value of the vacuum on the gauge
	b) Protection filter blocked or damaged	b) Replace the filter
	c) Connection tubes blocked, kinked or disconnected	c) Replace or reconnect the tubes, check the jar connections
	d) Shut-off valve blocked or damaged	d) Empty the jar, or disconnect the tube from the jar and unblock the shut-off valve. The unit will only work in the upright position
	e) Pump motor damaged	e) Refer to authorised service personnel
5. The float doesn't close	If the cap has been washed, ensure that the float is not partially detached	Insert the float into its place
6. The float doesn't close	The float it's covered by dirty material	Unscrewed the cap, leave the and put in on autoclave
7. Low suction	Foam inside the jar	Fill the jar to 1/3 full of ordinary water
<b>Faults 1 - 2 - 3 - 4 - 5 - 6 - 7</b>	<b>None of the remedies has achieved the desired results</b>	<b>Contact the seller or GIMA After-sales Assistance Service</b>

If the overflow security system it's activated, don't proceed with the liquid aspiration.

If the overflow security system doesn't work there are two cases:

1° case – If the overflow security system doesn't work the aspiration will be stopped by the bacteriological filter who avoid the liquid penetration inside the device.

2° case – If both the security system doesn't work, there is the possibility that liquid comes inside the device, in this case return the device to GIMA SpA technical service.

**The manufacturer will provide upon request electric diagrams, components list, descriptions, setting instructions and any other information that can help the technical assistance staff for product repair.**

**BEFORE EVERY CHECKING OPERATION, IN CASE OF ANOMALIES OR BAD FUNCTIONING, PLEASE CONTACT GIMA TECHNICAL SERVICE.  
GIMA SpA DOES NOT GIVE GUARANTEE IF INSTRUMENT, AFTER THE TECHNICAL SERVICE CHECKING, APPEARS TO BE TAMPERED.**

### **IMPORTANT INFORMATION FOR CORRECT DISPOSAL OF THE PRODUCT IN ACCORDANCE WITH EC DIRECTIVE 2002/96/EC:**

*In respect of art. 13 Decreto Legislativo 25 Luglio 2005, n. 151 "Attuation of European directives 2002/95/EC, 2002/96/EC and 2003/108/EC, for reduction in use of dangerous substances in the electric and electronic device and for garbage disposal"*

The symbol as over applied on the device or its packaging means that at the end of its useful life the product must not be disposed of with domestic waste. At the end of device useful, the user will must deliver it to the able collecting centres for electric and electronic garbage, or give back to the retailer in the moment of equivalent new device purchasing, one against one. Disposing of the product separately prevents possible negative consequences for the environment and for health, deriving from inadequate disposal. It also allows the recovery of materials of witch it's made up in order to obtain an important saving of energy and resources and to avoid negative effects to the ambient and health. In case of abusive disposal of device by user, will be applied administrative endorsements in compliance with current standard.


## INSTRUCTIONS

- Connect the short silicon tube **1**, with antibacterial filter **2**, to the suction connector **3**. The other tube, with one end connected to the filter must be connected with the other end to the jar's connector **4** where has been fixed the red float (security float).
- When the 90% of the volume of the jar is reached there is the activation of the security float (the float close the aspiration connector on the jar) to avoid liquid penetration inside the device. The device must be used on a plan of horizontal operation.
- Connect the long silicon tube **5** to the other jar's lid connector **6**
- Connect the other end of the long silicon tube **5** to the probe plastic connector **7**
- Connect the power cord to the device than connect the plug to the electrical mains supply.
- Unscrew the jar's lid and fill the jar 1/3 full of ordinary water (this for an easy cleaning operations and an rapid reaching of the functionality vacuum) then rescrew the lid on the jar correctly.
- Push switch **8** on position **I** to start suction
- Once finished push switch on **0** position and unplug
- To extract the accessories and start with cleaning it.

**WARNING:** Ensure that the IN marker on the filter is on the side facing the collection jar lid and fitted into the "VACUUM".  
A wrong connection causes immediate destruction in case of contact with sucked liquids.

### Filter assembling

The dual use version is available on request. Continuous use can be enabled by activating the illuminated switch. Alternatively, the pedal control can be used for aspirating on command, while leaving the operator's hands free. If you're using the pedal version, insert the plug of the pedal power cord into the appropriate socket located next to the equipment, push the second switch on the position **II**, so start and stop the suction operation by pressing relasind the pedal.

**Using FLOVAC® disposable collection system:** Before connecting the disposable collection system, remove the white ring fitted on the tank holder for a more comfortable insertion of the same container.

Connect the short tube with the throttle connector marked as "VACUUM" fitted on the lid of the disposable pocket. The container should be positioned on ON. The long silicon tube should be connected with the connector marked as "PATIENT".

In order to perform suction, close the connector marked as "TANDEM" with the lid provided.

If using the equipment with FLOVAC® disposable collection system antibacterial filters are not required since each disposable pocket is already provided with an inside filter. The hydrophobic antireflux and antibacterial filter provided will protect the equipment, while also acting as float valve deactivating vacuum generation when the maximum volume level is reached.

When the float valve intervenes signalling the device is too full, the suction source must be disconnected within no more than 5 minutes. **Warning:** **The accidental inversion of connections may cause contamination for the operator and/or for the vacuum generation equipment.**

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