GIMA

PROFESSIONAL MEDICAL PRODUCTS

Gima S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) Italy

Italia: tel. 199 400 401 - fax 199 400 403

Export: tel. +39 02 953854209/221/225 fax +39 02 95380056

gima@gimaitaly.com - export@gimaitaly.com

www.gimaitaly.com

SATURIMETRO OXY-6

OXIMETER OXY-6

OXYMÈTRE OXY-6

PULSOXIMETER OXY-6

SATURÓMETRO OXY-6

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Instructions to User

Read these instructions carefully before using this equipment. These instructions describe the operating procedures to be followed strictly. Failure to follow these instructions can cause measuring abnormality, equipment damage and personal injury. The manufacturer is NOT responsible for the safety, reliability and performance issues and any monitoring abnormality, personal injury and equipment damage due to user's negligence of the operation instructions. The manufacturer's warranty service does not cover such faults.

- The uncomfortable or painful feeling may appear if using the device ceaselessly, especially for the microcirculation barrier patients. It is recommended that the sensor should not be applied to the same finger for over 2 hours.
- For the individual patients, there should be a more prudent inspecting in the placing process. The device can not be clipped on the edema and tender tissue.
- The light (the infrared is invisible) emitted from the device is harmful to the eyes, so the user and the maintenance man, can not stare at the light.
- Testee can not use enamel or other makeup on the finger.
- Testee's fingernail can not be too long.
- Please peruse the relative content about the clinical restrictions and caution.

1. Safety

1.1 Instructions for Safe Operations

- Check the main unit and all accessories periodically to make sure that there is no visible damage that may affect patient's safety and monitoring performance. It is recommended that the device should be inspected once a week at least. When there is obvious damage, stop using the oximeter.
- Necessary maintenance must be performed by qualified service engineers ONLY. Users are not permitted to maintain it by themselves.
- The oximeter cannot be used together with the devices not specified in User's Manual.

1.2 Warnings

- *Explosive hazard—DO NOT use the oximeter in environment with inflammable gas such as some ignitable anesthetic agents.*
- *DO NOT use the oximeter while the testee measured by MRI and CT*
- *To dispose the device, the local law must be followed.*

1.3 Attentions

- *Keep the oximeter away from dust, vibration, corrosive substances, explosive materials, high temperature and moisture.*
- *If the oximeter gets wet, please stop operating it. When it is carried from cold environment to warm and humid environment, please do not use it immediately.*
- *DO NOT press the keys on front panel with sharp materials.*
- *High temperature or high pressure steam disinfection to the oximeter is not permitted. Refer to User's Manual for instructions of cleaning and disinfection.*
- *DO NOT have the oximeter immersed in liquid. When it needs cleaning, please wipe its surface with disinfect solution by soft material. Do not spray any liquid on the device directly.*

2. Overview

The pulse oxygen saturation is the percentage of HbO₂ in the total Hb of the blood, so-called the O₂ concentration in the blood. It is an important bio-parameter to the respiration. Many of the respiration disease will cause hypoxemia, even endanger the patient's health. As a result, monitoring the SpO₂ is indispensable in the clinical rescuing. The traditional method to measure SpO₂ is to analyze the sample of the patient's blood, so can get the partial pressure of oxygen and calculate the SpO₂ by use the blood-gas analyzer. This method is inconvenient and can not be used to monitor continuously. For the purpose of measuring the SpO₂ more easily

and accurately, GIMA developed the Fingertip Oximeter. The device can measure the pulse rate and blood perfusion index simultaneously. The Fingertip Oximeter is compact, convenient to use and carry and with low power consumption. You just need to put the fingertip into the sensor of the device, the SpO₂ value will appear on the screen immediately.

2.1 Features

- Small in size and lightweight;
- Color OLED, various display modes, display directions adjustable;
- The device can accurately measure SpO₂ value, pulse rate value and perfusion index;
- The device will automatically start measuring after putting finger into sensor;
- The device will power off automatically without signal for about 8 seconds;
- Audible & visual alarm
- Low voltage indication;

2.2 Major Applications and Scope

This device is applicable to home, hospital (including internal medicine, surgery, anesthesia, pediatrics, emergency room etc.), oxygen bar, the community medical center, alpine area and it also can be used before or after sports, and the like. The Fingertip Oximeter can detect SpO₂, pulse rate and blood perfusion index.

This device is not appropriate to be used for continuous monitoring.

2.3 Environment Requirements

Operating Temperature: 5°C ~40°C

Operating Humidity: 30%~80%

Atmospheric pressure: 70kPa~106kPa


2.4 Principle of Measurement

The measurement of pulse oximeter is that it uses a multi-functional oxyhemoglobinometer to transmit some narrow spectrum light bands through blood samples and to measure attenuation of spectrum with different wavelengths according to the characteristic that RHb, O2Hb, Met Hb and COHb absorb the light of different wavelength, thereby determining O2Hb saturation of different fractions. O2Hb saturation is called "fractional" O2Hb saturation.

$$\text{Fractional O2Hb saturation} = \left[\frac{\text{O2Hb}}{\text{RHb} + \text{O2Hb} + \text{Met Hb} + \text{COHb}} \right] \times 100$$

Oppositely, pulse oxygen oximeter measures functional O2Hb saturation:
 Functional O2Hb saturation = $\left[\frac{\text{O2Hb}}{\text{RHb} + \text{O2Hb}} \right] \times 100$

Present SpO2 oximeter transmits light of two wavelengths only, red light (wavelength 660 nm) and infrared (wavelength 940 nm), to differentiate HbO2 from HbR. One side of the sensor contains two LEDs, and the other side contains a photoelectric detector. SpO2 oximeter measures HbO2 saturation in the blood by the light plethysmograph when the pulse beats. The result is quite precise when HbO2 saturation is over 70% ~ 95%.

2.5 Caution

A. The finger should be placed properly (see the attached illustration of this manual), or else it may cause inaccurate measurement.

- B. The SpO₂ sensor and photoelectric receiving tube should be arranged in a way with the subject's arteriole in a position there between.
- C. The SpO₂ sensor should not be used at a location or limb tied with arterial canal or blood pressure cuff or receiving intravenous injection.
- D. Do not fix the SpO₂ sensor with adhesive or else it may result in venous pulsation and inaccurate measure of SpO₂.
- E. Make sure the optical path is free from any optical obstacles like rubberized fabric.
- F. Excessive ambient light may affect the measuring result. It includes fluorescent lamp, dual ruby light, infrared heater, direct sunlight and etc.
- G. Strenuous action of the subject or extreme electrosurgical interference may also affect the accuracy.
- H. Testee can not use enamel or other makeup on the finger.

3. Technical Specifications

A. Display mode: OLED Display

B. Power supply requirement:

1.5V (AAA size) alkaline battery x 2

Battery voltage: 2.7V~3.3V

C. Operating current: <40mA

D. SpO₂ Parameter Specifications:

Measuring range: 35%~99%

Accuracy: $\pm 2\%$ (during 75% - 99%)

Accuracy: $\pm 3\%$ (during 50% - 75%)

E. Pulse Rate Parameter Specifications:

Measuring range: 30bpm - 240bpm

Accuracy: ± 2 bpm or $\pm 2\%$ (which ever is greater)

F. Blood Perfusion Parameter Specifications:

Measuring range: 0.2%~20%

Accuracy: $\pm 0.1\%$ (during 0.2% - 2%)

Accuracy: $\pm 1\%$ (during 2% - 10%)

G. Resistance to interference of surrounding light:

The difference between the value measured in the condition of indoor natural light and that of darkroom is less than $\pm 1\%$.

H. Resistance to interference of man-made light:

Values of SpO₂ and Pulse Rate can be accurately measured by pulse oxygen simulator.

I. Dimensions: 66 mm (L) x 36 mm (W) x 33 mm (H)

Net Weight: 60g (including batteries)

J. Classification:

The type of protection against electric shock: Internally powered equipment.

The degree of protection against electric shock: Type BF equipment.

The degree of protection against harmful ingress of water: Ordinary equipment without protection against ingress of water.

Electro-Magnetic Compatibility: Group I, Class B

4. Accessories

A. A hanging cord

B. Two batteries

C. A pouch



D. A User Manual

E. Qualifies Certificate

5. Installation

5.1. Front view (figure 1)

5.2. Rear view (figure 2)


Figure 3

5.3. Battery

Refer to the figure, and insert the two AAA size batteries properly in the right direction.

Replace the cover.

- Please take care when you insert the batteries for the improper insertion may damage the device.

6. Operating Guide

1. Open the clip as shown in Figure 4.




Figure 4

2. Put finger onto the rubber cushion of the clip (make sure the finger is in the right position), and then clip the finger.

3. The device will power on automatically in 3 seconds, and start to display software version number.



4. Next enter into data display screen (as shown in Figure 5). After the waveform and value become stable, get the information directly from display screen.

*Figure 5**Figure 6*

5. When the display is shown as Figure 5, press Operation Key to switch display screen:

- Press Operation Key once, revolve the display screen 180°.
- Press Operation Key twice, the display screen will be changed as Figure 6.
- Press Operation Key three times, revolve the display screen (as shown Figure 6) 180°.
- Press Operation Key four times, the display screen will back to the screen as shown in Figure 5.

6. Longtime press Operation Key (about 2 seconds), the display is shown as Figure 7.

*Figure 7**Figure 8*

Difference between Figure 7 and Figure 5 are as follows:

- In Figure 5, SpO₂ and pulse rate are being monitored and displayed on the screen. In Figure 7, SpO₂ and poor perfusion index are being monitored and displayed on the screen.

7. When the screen displays as shown in Figure 7, press Operation Key circularly, the display screen will be switched between the two screens (as shown in Figure 7 and Figure 8), and four directions display alternately.

8. If without operation in 10 seconds on the screen which shows as Figure 7 or Figure 8, the display screen will return to the screen as shown in Figure 5 or Figure 6.

9. Alarm Indicator

When measuring, if SpO₂ value and pulse rate value exceeds the preset alarm limits, the device will alarm automatically and the value on the screen exceeding limit will blink; at this time press Display Key to shut down the alarm. Exceeding pulse rate alarm limit: sound twice as an interval;

Preset alarm range:

SpO₂ alarm: Low limit: 90%

Pulse Rate alarm: High limit 120bpm - Low limit: 50bpm

7. Repair and Maintenance

A. Please change the batteries when the low-voltage indicator lightens.

B. Please clean the surface of the device before using. Wipe the device with alcohol first, and then let it dry in air or wipe it dry.

C. Please take out the batteries if the oximeter will not be used for a long time.

D. The best storage environment of the device is -20°C to 55°C ambient temperature and less than 95% relative humidity.

*High-pressure sterilization cannot be used on the device.
Do not immerse the device in liquid.*

It is recommended that the device should be kept in a dry environment. Humidity may reduce the useful life of the device, or even damage it.

8. Troubleshooting

Trouble	Possible Reason	Solution
The SpO2 and Pulse Rate display instable	<ol style="list-style-type: none">1. The finger is not placed inside enough.2. The finger is shaking or the patient is moving.	<ol style="list-style-type: none">1. Place the finger properly and try again.2. Let the patient keep calm.
The device can not turn on	<ol style="list-style-type: none">1. The batteries are drained or almost drained.2. The batteries are not inserted properly.3. The device's malfunction	<ol style="list-style-type: none">1. Change batteries.2. Reinstall batteries.3. Please contact the local service center.
The indicator light is off suddenly	<ol style="list-style-type: none">1. The device will power off automatically when it gets no signal for 8 seconds.2. The batteries are almost drained.	<ol style="list-style-type: none">1. Normal.2. Change batteries.

Disposal: *The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment. For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.*

GIMA WARRANTY CONDITIONS

Congratulations for purchasing a GIMA product.

This product meets high qualitative standards both as regards the material and the production. The warranty is valid for 12 months from the date of supply of GIMA.

During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons. Labor costs and personnel traveling expenses and packaging not included.

All components subject to wear are not included in the warranty.

The repair or replacement performed during the warranty period shall not extend the warranty.

The warranty is void in the following cases: repairs performed by unauthorized personnel or with non-original spare parts, defects caused by negligence or incorrect use.

GIMA cannot be held responsible for malfunctioning on electronic devices or software due to outside agents such as: voltage changes, electro-magnetic fields, radio interferences, etc.

The warranty is void if the above regulations are not observed and if the serial code (if available) has been removed, cancelled or changed.

The defected products must be returned only to the dealer the product was purchased from. Products sent to GIMA will be rejected.

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